

**waterpik®  
waterflosser**

# Congratulations!

You are now on your way  
to better oral health!

## ✓ Register and Save

Register your product at [waterpik.com/registration](http://waterpik.com/registration) to receive over **\$50** in savings (includes **50%** off your next order of assorted tips + tip storage case).

## ✓ Write a Review

Creating innovative Waterpik® products that inspire better oral health is our passion. We invite you to leave a review by sharing your thoughts on the Cordless Express Water Flosser on [waterpik.com](http://waterpik.com) or the retailer's website.

## ✓ We're Here to Help

We hope you are loving your new Cordless Express Water Flosser. If for any reason you are not satisfied or have any product questions, please contact our fast and friendly Customer Service Team based in Fort Collins, Colorado. The Cordless Express Water Flosser is backed by a 1-year manufacturer's warranty.

Call: 1-800-525-2774  
Monday-Friday  
7 am - 5 pm MST

Chat or E-mail  
[service@waterpik.com](mailto:service@waterpik.com)



**Quick Start  
Guide Inside!**



## Register and Save!

### Save over \$50!

Register at [Waterpik.com/Registration](http://Waterpik.com/Registration) to receive exclusive discount code.

Offer includes  
**50% off**  
your next order  
of assorted  
tips + tip  
storage case.



### We'd Love to Hear from You!



Connect with Us  
**@Waterpik**



Call: 1-800-525-2774  
[www.waterpik.com](http://www.waterpik.com)

### Waterpik® 14 Day Challenge App




Make water flossing a habit! Try our new mobile App for the best water flossing experience.



# Quick Start Guide

## GETTING STARTED

- 1 Remove reservoir




- 2 To remove battery door, twist, and pull.



- 3 Insert 3 AA batteries.



- 4 To replace battery door, line-up arrows, push and twist.



- 5 Insert tip firmly until it clicks into place.



- 6 If the reservoir was removed, insert reservoir back on the unit. Fill reservoir with warm water.



## WATER FLOSS

- 7 With water switch in OFF position, place tip in mouth. Slide power switch to low (I) and increase to your liking.




- 8 Lean over sink and close lips enough to prevent splashing, while still allowing water to flow from mouth into the sink.



Keep unit upright during use to ensure continuous water flow.

## TECHNIQUE

- 9 Aim water at the gumline at a 90 degree angle. Follow the gumline and pause briefly between teeth.



**TO AVOID MESS: KEEP TIP IN MOUTH WHILE UNIT IS ON.**

## NEED HELP?

Visit [Waterpik.com/WelcomeWF](http://Waterpik.com/WelcomeWF)  
To watch a How To Video and find additional product support.

**Call Us! 1-800-525-2774**