

\$5 Rebate on Select Waterpik® Products

Offer expires December 31, 2025.

Valid on all Waterpik[®] Water Flosser and Complete Care Models Over \$45 retail.

Rebate not valid on individual Power Toothbrushes or Power Flossers. To receive your rebate check, please print your name, address, phone number and email address.

Go Online to waterpik.com/rebate to Submit your Rebate or Mail in the Following Form.

1. Completed Form

Name	
Home Address	
City	
State	Zip
Home Phone ()	

- 2. Original Sales Receipt. Circle the Waterpik® Product Purchased. (Must Mail with Form)
- 3. Cut out UPC (purchase seal) from Box.

(Must Mail with Form)



4. Mail the Collected Pieces to:

Inmar Rebate Center Attn: Patient Brochure \$5 MIR Offer Code WPIK75 PO BOX 426008 Laredo, TX 78044-6008

Items must be received by January 31, 2026.

Limit four rebates per household. Cannot be used in combination with other offers. Offer valid from 1/1/20-12/31/25. Valid for retail or professional purchases in the U.S. and possessions. Void where prohibited by law. Please allow 8–10 weeks for delivery. For inquiries, please call 1-877-322-8355 or go to waterpik.com/rebate.



CLINICALLY PROVEN RESULTS!

UP TO

EFFECTIVE PLAQUE REMOVAL

Number of Plaque from treated areas.

Treated

Treated

UP TO		
AS		
EFFECTIVE		

MORE EFFECTIVE THAN STRING FLOSS

Up to **2X** as effective for improving gum health vs. brushing and flossing.*





ESSENTIAL FOR IMPLANTS

Up to **2X** as effective for improving gum health around implants vs.string floss.*



Utilizina Plaaue Seeker™ tip

*Independent clinical study. Go to waterpik.com for details.

WHICH TIP IS RIGHT FOR YOU?

Classic Jet Tip	- General use	The Easy and Most Effective Way to Floss!™
Orthodontic Tip	- Braces - General use	Specially designed to reach underneath wires and around brackets to remove plaque.
Plaque Seeker™ Tip	- Implants - Bridges - Retainers - General use	Removes bacteria and debris from hard to reach areas around dental restorations.
Pik Pocket™ Tip	- Periodontal pockets - Furcations	Gently delivers medicaments and rinses deep below the gumline.
Implant Denture	- Implant retained	Specially designed to reach underneath and around implant retained

www.waterpik.com
©2022 Water Pik, Inc.

bacteria and debris

waterpik waterflosser



HEALTHIER GUMS AND BRIGHTER TEETH

Brushing alone is not enough. Despite the fact that virtually everyone brushes their teeth, half of all Americans have gum disease.

Clinical research confirms that the Waterpik® Water Flosser is significantly more effective than string floss, interdental brushes and air floss for removing plaque and improving gum health. That's why it is the #1 choice among dental professionals.

If you have Implants, Crowns, Bridges, Braces, or other dental work, the Waterpik® Water Flosser is essential.

SCAN TO DOWNLOAD OUR APP FOR THE BEST FLOSSING EXPERIENCE



waterpik

The Easy and Most Effective way to Floss

REBATE INSIDE



Learn More
about the Waterpik®
Water Flosser
and Waterpik®
Sonic-Fusion®

One Step to a **Complete Clean**



SUBSCRIBE & SAVE!

Save 20% + Free Shipping on Brush Head Subscriptions.



Scan to order replacement brush heads delivered right to your door!

waterpik sonic-fusion®



BRUSH AND FLOSS AT THE SAME TIME

- One Step to a complete Clean
- 2x as effective as traditional brushing & flossing
- 3 modes: Brush, Floss, Brush and Floss




Available in Dental Offices and Through Online Retailers













WHITENING

ION

Available in **Retail Stores**



2.0







Dental Professionals — to order more brochures please call 1-800-525-2020 or go to waterpik.com/brochures.



Waterpik® Sonic-Fusion®

Offer expires December 31, 2025.

Valid on Waternik® Sonic-Eusion® over \$149.99

To receive your rebate check, please print your name, address, phone number and email address

Go Online to waterpik.com/rebate to Submit your Rebate or Mail in the Following Form.

1. Completed Form

Name	
Home Address	
City	
State	Zip
Home Phone ()	
Email Address	

- 2. Original Sales Receipt. Circle the Waterpik® Product Purchased. (Must Mail with Form)
- 3. Cut out UPC (purchase seal) from Box. (Must Mail with Form)



4. Mail the Collected Pieces to:

Inmar Rebate Center Attn: Patient Brochure \$15 MIR Offer code WPIK77 PO BOX 426008 Laredo, TX 78044-6008

Save \$15 on purchase of One (1) Waterpik® Sonic-Fusion® after rebate. Offer valid on purchases made between 1/1/2020-12/31/2025. Must spend \$149.99 or more on One (1) Waterpik® Sonic-Fusion® to qualify for rebate. To receive your rebate, submit your Original Register receipt with the Qualifying Purchase, Purchase Price, and Purchase Date circled. Submit original UPC from product package. Must be received by 1/31/2026. Limit four rebates per household. Cannot be used in combination with other offers. Valid for retail or professional purchases in the U.S. and possessions. Void where prohibited by law.

Please allow 8-10 weeks for delivery. For inquiries, please call 1-877-322-8355 or go to waterpik.com/rebate.

