



Why just brush when you can get Complete Care?

Brushing alone is not enough, even with a high-end Sonic Toothbrush. Waterpik® Complete Care 5.0 combines Water Flosser and Sonic Toothbrush technologies to create a perfect daily oral care routine.

Brushing and flossing is now possible with one convenient device that saves counter space and power outlets.

waterpik®

waterpik® waterflosser®



ULTRA

Advanced features and performance

COMPLETE CARE 5.0

More effective brushing and flossing in one device!



CORDLESS PLUS

Cordless convenience

FOR KIDS

Easy-to-use, customized model for kids ages 6+



To find out where to purchase each item, go to waterpik.com

Dental Professionals — to order more brochures, please call 1-800-525-2020.

Waterpik®, Waterpik® (stylized), Plaque Seeker®, Pik Pocket™, The Easy and More Effective Way to Floss™, and WaterFlosser® (stylized) are trademarks of Water Pik, Inc.

Sonicare® is not a trademark of Water Pik, Inc.

©2016 Water Pik, Inc.
www.waterpik.com

waterpik®

The easy and more effective way to floss!™



COUPON \$5 INSIDE

**waterpik®
waterflosser®**

waterpik® waterflosser®



Brushing is Not Enough.

Despite the fact that virtually everyone brushes their teeth, half of all Americans have gum disease. Clinical research confirms that the Waterpik® Water Flosser is significantly more effective than string floss, Sonicare® Air Floss and Air Floss Pro for reducing gingivitis and improving gum health.¹ That's why it is the #1 choice among dental professionals.

And if you have dental work like implants, crowns, or braces, the Waterpik® Water Flosser is essential.

¹Independent clinical studies. Go to www.waterpik.com for details.

The Waterpik® Water Flosser: Clinically Proven Results

Effective Plaque Removal

Removes up to **99.9%** of plaque from treated areas.



Laboratory study, tooth surface magnified by electron microscope

Healthier Gums

Up to **50%** more effective for improving gum health vs. string floss.¹



Better Than Air Floss¹


Proven 80% more effective than Air Floss and 50% more effective than Air Floss Pro for improving gum health.



¹Independent clinical studies. Go to www.waterpik.com for details.



Clinically Proven Results

The Waterpik® Water Flosser provides a unique combination of water pressure and pulsations to clean deep between teeth and below the gumline, removing harmful bacteria and debris that traditional brushing and flossing can't reach.



| WHICH TIP IS RIGHT FOR YOU? | | |
|--|---|--|
| Classic Jet Tip  | - General Use | The Easy and More Effective Way to Floss! [™] |
| Toothbrush Tip  | - General Use | The convenience of brushing and Water Flossing at the same time. |
| Plaque Seeker® Tip  | - Implants - Crowns - Bridges - General Use - Retainers | Removes bacteria and debris from hard to reach areas around dental restorations. |
| Orthodontic Tip  | - Braces - General Use | Specially designed to reach underneath wires and around brackets to remove plaque. |
| Pik Pocket™ Tip  | - Periodontal Pockets - Furcations | Gently delivers medicaments and rinses deep below the gumline. |

5 UNIQUE TIPS AVAILABLE

\$5 Mail-in Rebate on Select Waterpik® Water Flossers

Offer expires December 31, 2018

Valid only on the purchase of the Waterpik® Ultra, Cordless Plus, Complete Care and For Kids models.

Rebate not valid on Power Toothbrushes or Power Flossers. To receive your rebate check, please print your name, address, phone number and email address.

Name _____

Home Address _____

City _____

State _____ Zip _____

Home Phone (_____) _____

Email Address _____

Mail the original, completed rebate offer form, original sales receipt with Waterpik® Product purchase circled, and original proof of purchase seal (UPC). Each submission must include 1 form to qualify.

Inmar Rebate Center
Attn: Patient Brochure \$5 MIR
Offer Code WPIK36
P.O. Box 426013
Del Rio, TX 78842-6013

Items must be received by January 31, 2019

Limit four rebates per household. Cannot be used in combination with other offers. Offer valid from 4/1/16 to 12/31/18. Valid for retail or professional purchases in the U.S. and possessions. Void where prohibited by law. Please allow 8-10 weeks for delivery. For inquiries, please call 1-877-322-8355 or go to www.inmarrebates.com.

FN 20021840-SSH-F AA

waterpik®



Results Guaranteed

Healthier gums in just 14 days guaranteed or your money back.

If after 14 days of daily use you are not completely satisfied that this product has reduced the symptoms associated with gingivitis (e.g., inflammation and bleeding of gums) return it for a full refund, no questions asked.