






WATERPIK® WATER FLOSSER TIP GUIDE

WHICH TIP IS RIGHT FOR YOU?

<input type="checkbox"/> Classic Jet Tip 	- General use	The Easy and Most Effective Way to Floss!™
<input type="checkbox"/> Orthodontic Tip 	- Braces - General use	Specially designed to reach underneath wires and around brackets to remove plaque.
<input type="checkbox"/> Plaque Seeker™ Tip 	- Implants - Bridges - Retainers - General use	Removes bacteria and debris from hard to reach areas around dental restorations.
<input type="checkbox"/> Pik Pocket™ Tip 	- Periodontal pockets - Furcations	Gently delivers medicaments and rinses deep below the gumline.
<input type="checkbox"/> Implant Denture Tip 	- Implant retained dentures	Specially designed to reach underneath and around implant retained dentures to remove bacteria and debris.

Instructions: