

If it doesn't say **waterpik™**... it isn't **waterpik™**



FOR SAFETY
AND EFFICACY



WHICH PRODUCT IS RIGHT FOR YOU?

COUNTERTOP MODELS

- Excellent Performance
- Adjustable water pressure settings
- Includes multiple Tips for different patients needs

Great for patients with plaque and gingivitis

TOOTHBRUSH + WATER FLOSSER

- Sonic-Fusion offers convenience of a toothbrush with a built-in water flosser
- Senonic Complete Care offers a separate toothbrush and water flosser for a customized experience

Great for those with limited counterspace and 1 outlet use

HANDHELD MODELS

- Ultimate Convenience
- Lightweight and portable
- Waterproof, can be used in the shower

Great for travelers and caregivers

WATERPIK™ REBATES

Offer expires December 31, 2025



Scan or go to
waterpik.com/rebate
to submit your rebate

Waterpik™ Sonic-Fusion™ or Senonic™ Complete Care

\$15^{USD} Rebate WPIK77

Valid on product over \$149.99

Waterpik™ Select Products

\$5^{USD} Rebate WPIK75

Valid on water flossers over \$45 retail


Rebate only valid on select Waterpik water flosser retail valued over \$45 purchased at retail. To receive your rebate, submit your original receipt with the Qualifying Purchase, Purchase Price, and Purchase Date circled. Submit original UPC from product package. Limit 4 rebates per household. For redemption or for inquiries, please call 1-877-322-8355 or go to waterpik.com/rebate. Each submission must include 1 form per individual product. Cannot be used in combination with other offers. Offer good until 12/31/25. Items must be received by January 31, 2026. Valid on retail purchases in USD. Void where prohibited by law. Please allow 8-10 weeks for delivery. Unless expressly prohibited by law, payee authorizes reasonable dormancy fees deducted if check not cashed within 180 days. Reproduction of this form is expressly prohibited. Any other use of this form constitutes fraud.

UP TO
99.9%
PLAQUE
REMOVAL
from treated areas

EFFECTIVE PLAQUE REMOVAL

Removes up to **99.9%** of plaque from treated areas!¹

Laboratory study, tooth surface magnified by electron microscope



Healthier Gums and Cleaner Teeth!

Waterpik Delivers Clinically Proven Results!

MORE EFFECTIVE THAN STRING FLOSS

Up to **2X** as effective for improving gum health vs. brushing and flossing!¹

¹Independent clinical study. Go to waterpik.com for details.

UP TO
2X
AS
EFFECTIVE

WHICH TIP IS RIGHT FOR YOU?

	<input type="checkbox"/> Precision Tip		- General use	The Easy and Most Effective Way to Floss! TM <small>*compared to string floss and interdental brushes</small>
	<input type="checkbox"/> Orthodontic Tip		- Braces - General use	Specially designed to reach underneath wires and around brackets to remove plaque.
	<input type="checkbox"/> Plaque SeekerTM Tip		- Implants - Bridges - Retainers - General use	Removes bacteria and debris from hard to reach areas around dental restorations.
	<input type="checkbox"/> Pik PocketTM Tip		- Periodontal pockets - Furcations	Gently delivers medicaments directed by a dental professional and rinses deep below the gumline.
	<input type="checkbox"/> Implant Denture Tip		- Implant retained dentures	Specially designed to reach underneath and around implant retained dentures to remove bacteria and debris.



Need Help or New to Water Flossing.

Scan this app for the best experience at home