

# Precision oral care for healthier gums\* in just 2 weeks

## Which product is right for you?



### ☐ Cordless handheld models

- **Great for small spaces and caregivers**
- Rechargeable, portable design
- 360 tip rotation
- 3 water pressure settings
- May be used in the shower



### ☐ Countertop performance models

- **Recommended for those with counter space**
- Maneuverable handle with 360 tip rotation
- 10 water pressure settings
- 90+ seconds of flossing time on low setting



### ☐ Flossing toothbrush

- **Great for those with limited counterspace and outlets**
- Brush, water floss, or both with the touch of a button

\*compared to string floss

Rebates  
redemptions are  
available through  
these retailers:



## WATERPIK™ REBATES

Not compatible with the Waterpik™ Patient Discount Program



Scan or go to  
[waterpik.com/rebate](https://waterpik.com/rebate)  
to submit your rebate

☐ Waterpik™ Sonic-Fusion™ or Sensonic™ Complete Care  
**\$15 Rebate WPIK77** Valid on product over \$149.99

☐ Waterpik™ Select Products  
**\$5 Rebate WPIK75** Valid on water flossers over \$45 retail

Rebate only valid on any Waterpik water flosser retail valued over \$45 purchased in retail. To receive your rebate, submit your Original Register receipt with the Qualifying Purchase, Purchase Price, and Purchase Date circled. Submit original UPC from product package. Limit 2 rebates per household, name & address. Each submission must include 1 form per individual product. Cannot be used in combination with other offers. Offer good until 12/31/2028 and items must be received January 31, 2029. Valid on retail purchases in the US only. Void where prohibited by law. Please allow 8-10 weeks for delivery. Unless expressly prohibited by law, payee authorizes reasonable dormancy fees deducted if check not cashed within 180 days. Reproduction of this form is expressly prohibited. Any other use of this form constitutes fraud.

# Waterpik™ Delivers *Proven Results!*

## More Effective Than String Floss

Up to **2X** as effective for improving gum health vs. brushing and flossing.<sup>1</sup>

<sup>1</sup>Independent clinical study. Go to [waterpik.com](http://waterpik.com) for details.


UP TO  
**2X**  
AS  
EFFECTIVE

UP TO  
**99.9%**  
PLAQUE  
REMOVAL  
from treated  
areas

## Effective Plaque Removal

Removes up to **99.9%** of plaque from treated areas.<sup>1</sup>

Laboratory study, tooth surface magnified by electron microscope













***Your new water flosser may be eligible for FSA, HRA, or HSA reimbursement!***

Visit [waterpik.com/FSA](http://waterpik.com/FSA) to print a letter of medical necessity

## Which tip is right for you?

Included with

	<input type="checkbox"/> <b>Precision Tip</b> 	- General use	Designed for everyday use to clean deep between teeth and below the gumline.	Every Waterpik™ water flosser
	<input type="checkbox"/> <b>Orthodontic Tip</b> 	- Braces - General use	Specially designed to reach underneath wires and around brackets to remove plaque.	Aquarius™, Ultra, ION, Cordless Advanced
	<input type="checkbox"/> <b>Plaque Seeker™ Tip</b> 	- Implants - Bridges - Retainers - General use	Removes bacteria and debris from hard to reach areas around dental restorations.	Aquarius™, Ultra, ION, Cordless Advanced
	<input type="checkbox"/> <b>Pik Pocket™ Tip</b> 	- Periodontal pockets - Furcations	Gently delivers medicaments directed by a dental professional and rinses deep below the gumline.	Aquarius™, Ultra, ION
	<input type="checkbox"/> <b>Implant Denture Tip</b> 	- Implant retained dentures	Specially designed to reach underneath and around implant retained dentures to remove bacteria and debris.	ION



**New to water flossing or need help?**  
Scan for the best experience at home.