How to Use a Waterpik® Water Flosser



Make sure to use warm water in the reservoir to avoid sensitivity

Place the flossing tip in your mouth before the unit on



HELPFUL TIPS

Mouthwash such as TheraBreath can be used in water flossers. Flush with water after use.

HELPFUL TIPS

New Users: Start on lowest setting and pressure

Scan to download for a **Better Flossing Experience**



Which Tip is Right for Your Patient?

which tip is kight for four Patients							
Tip	Usage	Benefit	Recommend for use with:				
Classic Jet Tip JT-100E, JT-450E	- General use	The Easy and Most Effective Way to Floss!™	All Waterpik Water Flossers				
Orthodontic Tip OD-100E	- Braces - General use	Specially designed to reach underneath wires and around brackets to remove plaque.	All Waterpik Water Flossers				
Plaque Seeker™ Tip PS-100E	- Implants - Bridges - Retainers - General use	Removes bacteria and debris from hard to reach areas around dental restorations.	All Waterpik Water Flossers				
Pik Pocket™ Tip	- Periodontal pockets - Furcations	Gently delivers medicaments and rinses deep below the gumline.	Countertop or ION Water Flossers				
Tongue Cleaner Tip TC-100E	- Fresher breath	Removes bacteria and sulfur compounds that reside on the tongue and cause bad breath.	All Waterpik Water Flossers				
Implant Denture Tip DT-100E	- Implant retained dentures	Specially designed to reach underneath and around implant retained dentures to remove bacteria and debris.	Countertop or ION Water Flossers				

Where to buy a Waterpik Water Flosser or Accessories























waterpik	Which Waterpik® Water Flosser is right for me?						
Туре	Performance Water Flosser	Performance Water Flosser	Performance Water Flosser	Handheld Water Flosser	Handheld Water Flosser	Flossing Toothbrush	
ADA Accepted American Dental Association ®	Aquarius Professional	Ultra	ION Professional	Cordless Advanced Professional	Cordless Express	Sonic-Fusion® 2.0 Professional	
Best for	 Easily reaches hard to clean areas like back teeth Convenient 1-minute timer and 30-second pacer Multiple users 	 Easily reaches hard to clean areas like back teeth Multiple users 	 Full performance with cordless convenience Bathrooms with limited outlets 	Small bathroomsIn-shower useTravelRe-chargeable	 Small bathrooms In-shower use Travel Comes with 3 AA batteries 	 Ultimate in convenience Brush and floss with one device Global voltage 	
Reservoir Capacity	90 seconds; can do the entire mouth	90 seconds; can do the entire mouth	90 seconds; can do the entire mouth	45 seconds; fill twice for entire mouth	30 seconds; fill twice for entire mouth	60+ seconds; can do the entire mouth	
Pressure Settings	10	10	10	3	2	10 2 brush speeds	
Accessories Included	Tips Classic Jet Tip (JT-100E) Pik Pocket Tip (PP-100E) Plaque Seeker Tip (PS-100E) Orthodontic Tip (OD-100E) Toothbrush Tip (TB-100E) Tip storage case	6 Tips 2 Classic Jet Tip (JT-100E) 1 Pik Pocket Tip (PP-100E) 1 Plaque Seeker Tip (PS-100E) 1 Orthodontic Tip (OD-100E) 1 Toothbrush Tip (TB-100E)	• 7 Tips 2 Classic Jet Tip (JT-100E) 1 Pik Pocket Tip (PP-100E) 1 Plaque Seeker Tip (PS-100E) 1 Orthodontic Tip (OD-100E) 1 Implant Denture Tip (DT-100E) 1 Tongue Cleaner (TC-100E) • Tip storage case	4 Tips 2 Classic Jet Tip (JT-450E) 1 Plaque Seeker Tip (PS-100E) 1 Orthodontic Tip (OD-100E) Tip storage case	• 2 Tips 2 Classic Jet Tip (JT-450E)	 2 Patented flossing brush heads 1 Compact size with cover 1 Full-size with cover 	
Warranty	3-Year	3-Year	3-Year	2-Year	1-Year	3-Year	
Price	\$\$	\$\$	\$\$	\$\$	\$	\$\$\$	