





**waterpik®**

**ORAL HEALTH PRODUCT GUIDE**



# waterpik® sonic-fusion® 2.0

PROFESSIONAL



### CLINICALLY PROVEN

Removes up to **99.9%** of plaque from treated areas.<sup>1</sup>

Up to **2X** as effective for removing plaque vs. regular brushing and flossing.<sup>1</sup>

### WORLD'S 1ST FLOSSING TOOTHBRUSH

Waterpik® Sonic-Fusion® seamlessly switches from brushing to flossing, and removes the harmful bacteria and debris deep between teeth and below the gumline that traditional brushing and flossing can't reach.

- 3 Modes - Brush, Floss, Brush + Floss
- Clinically proven
- Up to 2X as effective as regular brushing and flossing
- Limited three-year warranty



Includes:  
1 Compact Flossing Brush Head  
1 Full Size Flossing Brush Head  
and toothbrush case



# waterpik® ION PROFESSIONAL waterflosser™



**CLINICALLY PROVEN**

Removes up to **99.9%** of plaque from treated areas.<sup>1</sup>

Up to **50%** more effective for improving gum health vs. string floss.



**COUNTERTOP PERFORMANCE. CORDLESS CONVENIENCE.**

- Same superior pressure as the Aquarius countertop
- Advanced pressure control system with 10 settings
- Reservoir capacity of 90 seconds
- Unit stays charged for up to 4 weeks
- Magnetic Handle Cradle
- Seven water flossing tips included
- Limited three-year warranty

7 TIPS INCLUDED	BEST FOR
Classic Jet Tip	- General use
Plaque Seeker™ Tip	- Implants - Retainers - Crowns - General use - Bridges
Orthodontic Tip	- Braces - General use
Pik Pocket™ Tip	- Periodontal pockets - Furcations
Tongue Cleaner Tip	- Fresher Breath
Implant Denture Tip	- Fixed implant bridges and dentures



# waterpik® waterflosser™

CORDLESS  
ADVANCED



**CLINICALLY  
PROVEN**




Removes up to **99.9%** of plaque from treated areas.<sup>1</sup>

Up to **50%** more effective for improving gum health vs. string floss.

**WATERPROOF HANDHELD DESIGN.  
PERFECT FOR THE SHOWER.**

Cordless Advanced features a magnetic 4-hour rapid charging system, ultra quiet operation, waterproof design, and global voltage. It's perfect for small spaces, travel, and use in the shower. Available in several designer colors.

- Portable and lightweight
- Perfect for travel and small bathrooms
- 45 Seconds of water capacity
- 4-hour rapid charge
- 3 pressure settings
- Extra-quiet design
- Tip storage case include
- Waterproof and may be used in the shower
- Global voltage compatible
- 360-degree tip rotation
- Microfiber travel bag and travel water plug included to prevent leakage during transport
- Four water flossing tips included
- Two-year warranty

4 WATER FLOSSER TIPS	
	<p><b>Classic Jet Tip</b></p> <p>- General use</p>
	<p><b>Plaque Seeker Tip</b></p> <p>- Implants - Crowns - Bridges - Retainers - General use</p>
	<p><b>Orthodontic Tip</b></p> <p>- Braces - General use</p>



# waterpik® waterflosser™

AQUARIUS®  
PROFESSIONAL








**CLINICALLY PROVEN**

Removes up to **99.9%** of plaque from treated areas.<sup>1</sup>

Up to **50%** more effective for improving gum health vs. string floss.

**THE EASY AND MOST EFFECTIVE WAY TO FLOSS™**

- 10 Pressure settings (0.7 to 7.0 bar/ 10-100 psi)
- Advanced pressure control system
- Reservoir capacity of 90+ seconds
- On/Off water control on handle
- 2 Modes - Floss Mode and Hydro-Pulse Massage Mode
- 1-minute timer with 30-second pacer
- LED information panel
- Seven water flossing tips included
- Limited three-year warranty

7 TIPS INCLUDED	BEST FOR
 <b>Classic Jet Tip</b>	- General use
 <b>Plaque Seeker™ Tip</b>	- Implants - Retainers - Crowns - General use - Bridges
 <b>Orthodontic Tip</b>	- Braces - General use
 <b>Pik Pocket™ Tip</b>	- Periodontal pockets - Furcations
 <b>Toothbrush Tip</b>	- General use

# waterpik® waterflosser®

WHITENING  
PROFESSIONAL



IN-HANDLE  
WHITENING  
INFUSER

Actual Results



Before:  
Accumulation  
of tea stains



4 weeks of  
Waterpik®  
Whitening

CLINICALLY  
PROVEN




Removes up to **99.9%** of plaque from treated areas.<sup>1</sup>

Up to **25%** stain reduction after 2 weeks.<sup>1</sup>

## NATURALLY WHITEN WHILE YOU FLOSS

The perfect product for naturally whiter teeth and fresh breath, with all the benefits of water flossing. It restores your natural whiteness in just 4 weeks. The Whitening infuser technology expertly mixes water with a gentle stain removal agent that safely restores your natural whiteness, removes plaque, and promotes healthier gums.

- 10 Pressure settings
- Reservoir capacity of 90+ seconds
- On/Off water control on handle
- In-Handle Whitening Infuser
- 2 Modes – Floss Mode and Hydro-Pulse Massage Mode
- Four whitening water flosser tips included
- Includes 30 Mint Fresh Whitening Tablets
- Limited three-year warranty

4 TIPS INCLUDED	BEST FOR
 <p>Classic Jet Tip</p>	- General use
 <p>Orthodontic Tip</p>	- Braces - General use
 <p>Plaque Seeker™ Tip</p>	- Implants - Crowns - Bridges - Retainers - General use

AVAILABLE AT WATERPIK.COM



## Dispense Waterpik® Products

Dispensing the Waterpik® Water Flosser is a great way to get an effective product in your patient's hands quickly and get them on their way to healthier gums and teeth! Whether you plan to sell the product or include it in a treatment plan, Waterpik® offers several resources to help you be successful.



### sonic-fusion® 2.0

#### PROFESSIONAL

- 2 minute brushing timer with 30 second pacer
- Advanced pressure control
- Includes two Flossing toothbrush heads
- Recharge indicator light
- Deluxe toothbrush travel case
- Limited three-year warranty
- \$30 Patient Mail-in Rebate



#### AQUARIUS® PROFESSIONAL

- Our most advanced Water Flosser
- On/Off water control on handle
- 2 Water Flossing Modes — Floss and Massage
- Includes 7 Water Flosser Tips — for all patient types
- 1 minute timer with 30-second pacer
- Limited three-year warranty
- \$20 Patient Mail-in Rebate



#### ION PROFESSIONAL

- Same superior performance as Aquarius counter top
- Advanced pressure control system with 10 settings
- Reservoir capacity of 90 seconds
- Unit stays charged for up to 4 weeks
- Magnetic Handle Cradle
- Seven water flossing tips included
- Limited three-year warranty
- \$20 Patient Mail-in Rebate



#### CORDLESS ADVANCED

- Re-chargeable — no batteries required
- 4 Hour rapid charge system
- LED pressure and charge indicator
- Waterproof for flossing in the shower
- 3 Pressure settings
- Includes 4 Water Flosser tips
- Travel Bag and Tip Storage Case
- Limited two-year warranty
- \$20 Patient Mail-in Rebate





**WHITENING PROFESSIONAL**



- In-Handle Whitening Infuser
- 90 Seconds of Water Capacity
- 30 Mint Fresh Whitening Tablets
- Includes 4 Water Flosser Tips
- Accepted by the American Dental Association
- Limited three-year warranty

<sup>1</sup>Independent clinical study. Go to [waterpik.com](http://waterpik.com) for details.

**waterpik<sup>®</sup>**  
**PRO CONNECT<sup>SM</sup>**

Waterpik<sup>®</sup> Pro Connect is the go-to source for ordering trial and dispensing products for your office. Sign-up at [pro.waterpik.com](http://pro.waterpik.com) to receive exclusive dental professional pricing on Waterpik<sup>®</sup> Water Flossers and Waterpik<sup>®</sup> Sonic-Fusion<sup>®</sup> and order Patient Education Brochures.



**Professional Recommending Guide**

Dispensing Waterpik<sup>®</sup> Water Flossers is a great way to get an effective product in your patient's hands quickly and get them on their way to healthier gums and teeth! Whether you plan to sell the product or include it in a treatment plan, Waterpik<sup>®</sup> offers several resources to help you be successful.

[View recommending guide on waterpik.com](http://waterpik.com)



# waterpik®

## Professional series

Available at dental offices and online



Sonic Fusion® 2.0



Aquarius®  
Professional



ION  
Professional



Whitening  
Professional



Cordless  
Advanced

# waterpik®

## Retail series

Available at retail stores



Sonic Fusion® 2.0



Ultra



ION



Complete Care 5.0  
Sonic



Whitening



Cordless Plus



Cordless Express



Nano

# Health Benefits of Water Flossing

The Waterpik® Water Flosser is clinically proven to improve gum health and remove plaque and bacteria from hard to reach areas.



## EFFECTIVE PLAQUE REMOVAL

Removes up to **99.9%** of plaque from treated areas.<sup>1</sup>




## HEALTHIER GUMS

Reduced bleeding by up to 93% in 4 weeks.<sup>1</sup>



## PROVEN SAFE

The Waterpik® Water Flosser has been **proven safe** based on over 50 years of clinical studies.<sup>1</sup>



## SUPERIOR CLEANING AROUND BRACES

Up to **3X** as effective for removing plaque around braces vs. string floss.<sup>1</sup>



## ESSENTIAL FOR IMPLANTS

Up to **2x** as effective for improving gum health around implants vs. string floss.<sup>1</sup>



## REDUCES SUBGINGIVAL BACTERIA

Improves microflora and reduces counts of pathogenic bacteria in pockets up to 6 mm in depth.<sup>1</sup>



## IMPLANT DENTURE TIP

Specially designed to reach underneath and around implant retained dentures to remove bacteria and debris.

<sup>1</sup>Independent clinical study.  
Go to [waterpik.com](http://waterpik.com) for details.

# Clinical Research Highlights




All Water Pik, Inc. clinical research is conducted by independent universities and research facilities and is published in peer reviewed journals. Recent research proves conclusively that the Waterpik® Water Flosser is more effective than string floss for removing plaque and reducing gingivitis and gingival bleeding.

## SONIC-FUSION®

### The Waterpik® Sonic-Fusion®: Up to 2X as Effective as Regular Brushing & Flossing for improving Gingival Health

#### Comparison of sonic-flossing toothbrush to brushing and flossing on inflammation

Lyle DM, Qaqish JG, Goyal CR, Schuller R.  
IADR/AADR/CADR General Session & Exhibition.  
Presented at IADR Virtual Oral Session, July 24, 2021.  
Submitted for publication.






\*Statistically significant difference, p<0.001

## SONIC-FUSION®

# Waterpik® Sonic-Fusion®: Up to Twice as Effective as Regular Brushing and Flossing for Removing Plaque and Improving Gingival Health.

Comparison of a novel sonic toothbrush to standard brushing and flossing on clinical signs of inflammation: A randomized clinical trial

Goyal CR, Qaqish, JG, Schuller R, Lyle D. Compend Contin Educ Dent. 2018;39(suppl 2):14-22.



\*whole mouth results.  
Statistically significant, p<0.001


## SONIC-FUSION®

### Waterpik® Sonic-Fusion®: Significantly More Effective than Sonicare® DiamondClean for Removing Plaque and Improving Gingival Health


Direct comparison of a novel sonic toothbrush to a traditional sonic brush on clinical signs of inflammation: A randomized clinical trial

Goyal CR, Qaqish, JG, Schuller R, Lyle D. Compend Contin Educ Dent. 2018;39(suppl 2):14-22.


Gingival Bleeding Reduction\*



Gingival Inflammation Reduction\*



Plaque Reduction\*



\*whole mouth results.  
Statistically significant, p<0.001


## WATER FLOSSER

# The Waterpik® Water Flosser Significantly Improves Oral Health Benefits over a High-End Oscillating Electric Toothbrush

### Efficacy of the Use of a Water Flosser in Addition to an Electric Toothbrush on Clinical Signs of Inflammation: 4-Week Randomized Controlled Trial


Goyal CR, Qaqish JG, Schuller R, Lyle DM. Journal 2019; Compend Contin Ed Dent 2020, 41(3): March 2020, Volume 41, Issue 1.

Gingival Bleeding Reduction\*




P<0.001

Gingival Inflammation Reduction\*



P<0.001

Plaque Reduction\*



p=0.003

\*whole mouth results


Statistically significant difference

## WHITENING

### Waterpik® Whitening Water Flosser: Improved Stain Removal Over Tooth Brushing Alone

Evaluation of Tooth Whitening using a Liquid Dentifrice Delivered by the Waterpik® Whitening Water Flosser.

Milliman JL, Milleman KR, Battershell K, Lyle DM. Study conducted at Salus Research, Fort Wayne, IN. 2014.



Before.


After 2 Weeks of Using the Waterpik® Water Flosser with a Liquid Dentifrice.

## WATER FLOSSER

### Waterpik® Water Flosser Removes 99.9% of Plaque Biofilm After 3-Second Treatment

Biofilm Removal with a Dental Water Jet

Gorur A, Lyle DM, Schaudinn C, Costerton JW. Compend Contin Ed Dent 2009; 30 (Suppl 1):1-6. Study conducted at the University of Southern California School of Dentistry, USC Center for Biofilms, Los Angeles, California.






## WATER FLOSSER

### Waterpik® Water Flosser: Significantly More Effective than String Floss for Removing Plaque

Evaluation of the Plaque Removal Efficacy of a Water Flosser Compared to String Floss in Adults After a Single Use

Goyal CR, Lyle DM, Qaqish JG, Schuller R. J Clin Dent 2013; 24(2):37-42. Study conducted at BioSci Research Canada, Ltd., 3.




\*Statistically significant difference, p<0.001

### Waterpik® Water Flosser: Up to Twice as Effective as String Floss for Reducing Gingival Bleeding

The Effect of Different Interdental Cleaning Devices on Gingival Bleeding

Rosema NAM, et al. J Int Acad Periodontol 2011; 13(1):2-10. Study conducted at the University of Amsterdam, Academic Center for Dentistry, Amsterdam.




\*Statistically significant difference, p<0.05

## WATER FLOSSER

### Waterpik® Water Flosser: Up to 50% More Effective than String Floss for Reducing Gingivitis

Comparison of Irrigation to Floss as an Adjunct to Toothbrushing: Effect on Bleeding, Gingivitis and Supragingival Plaque

Barnes CM, Russell CM, Reinhardt RA et al. J Clin Dent, 2005; 16(3): 71-77. Study conducted at the University of Nebraska Medical Center, College of Dentistry, Lincoln, Nebraska.




\*Statistically significant difference, p=0.045

### The Waterpik® Water Flosser: Significantly More Effective than Interdental Brushes for Improving Gingival Health

Comparison of Water Flosser and Interdental Brush on Reduction of Gingival Bleeding and Plaque: A Randomized Controlled Pilot Study.

Goyal CR, Lyle DM, Qaqish JG, Schuller R. J Clin Dent 2016; 27: 61-65.



Whole Mouth


\*Statistically significant difference, p<0.001

## WATER FLOSSER

### Waterpik® Water Flosser: Up to 3 Times as Effective as String Floss for Orthodontic Patients

The Effect of a Dental Water Jet with Orthodontic Tip on Plaque and Bleeding in Adolescent Orthodontic Patients with Fixed Orthodontic Appliances


Sharma NC, Lyle DM, Qaqish JG et al. Am J Orthod Dentofacial Orthop 2008; 133(4): 565-571. Study conducted at BioSci Research Canada, Ltd., Mississauga, Ontario, Canada.









### Waterpik® Water Flosser: Up to 2 Times as Effective as String Floss for Implant Patients

Comparison of the Effect of Two Interdental Cleaning Devices Around Implants on the Reduction of Bleeding: A 30-day Randomized Clinical Trial

Magnuson B, Harsono M, Stark PC, et al. Compend Contin Ed Dent 2013; 34(Special Issue 8):2-7. Study conducted at Tufts University, School of Dental Medicine, Boston, Massachusetts.



## WATER FLOSSER TIP SELECTOR GUIDE

WHICH TIP IS RIGHT FOR YOUR PATIENT?		
<b>Classic Jet Tip</b> 	<ul style="list-style-type: none"> <li>- General use</li> </ul>	The Easy and Most Effective Way to Floss!™
<b>Orthodontic Tip</b> 	<ul style="list-style-type: none"> <li>- Braces</li> <li>- General use</li> </ul>	Specially designed to reach underneath wires and around brackets to remove plaque.
<b>Plaque Seeker™ Tip</b> 	<ul style="list-style-type: none"> <li>- Implants</li> <li>- Bridges</li> <li>- Retainers</li> <li>- General use</li> </ul>	Removes bacteria and debris from hard to reach areas around dental restorations.
<b>Pik Pocket™ Tip</b> 	<ul style="list-style-type: none"> <li>- Periodontal pockets</li> <li>- Furcations</li> </ul>	Gently delivers medicaments and rinses deep below the gumline.
<b>Tongue Cleaner Tip</b> 	<ul style="list-style-type: none"> <li>- Fresher Breath</li> </ul>	Removes bacteria and sulfur compounds that reside on the tongue and cause bad breath.
<b>Implant Denture Tip</b> 	<ul style="list-style-type: none"> <li>- Implant retained dentures</li> </ul>	Specially designed to reach underneath and around implant retained dentures to remove bacteria and debris

### Tip Storage Case

This convenient storage case holds 6 Waterpik® Water Flosser Jet Tips — it includes our most popular general use tips.

**Full Case Includes:**

- 2 Classic Jet Tips
- 2 Plaque Seeker™ Tips
- 2 Orthodontic Tips



## SUPPORT MATERIALS

[Click Here for Office Resources](#)

### Water Pik, Inc. — Supporting Your Patients and Improving the Bottom Line of Your Practice

We provide complimentary Continuing Education courses for you and your staff, and free educational materials for your patients. Our professional pricing and quantity discount programs make dispensing the Waterpik® Water Flosser a highly profitable opportunity for any dental practice.




#### Complimentary Practice Support Materials:

- Patient Education Brochures and Holders
- Professional and Patient Education Videos (online)
- Clinical Research Abstracts

#### Continuing Education Self-Study Courses with credits (Available at [waterpik.com](http://waterpik.com))

- **The Water Flosser:** Your Tool For Optimal Oral Health
- **Sonic-Fusion®:** Changing the Game in Interdental Care.
- **Risk-Based Periodontal Therapy:** A 21st-Century Perspective
- **Flossing Unraveled:** The Evidence on Interdental Cleaning
- **Vaping, Nicotine, and Tobacco:** New Trends, New Concerns



#### Online Dental Professional Newsletter —

- Sign up at [professional.waterpik.com](http://professional.waterpik.com)


#### Professional Pricing and Quantity Discounts

- Contact us at **1-800-525-2020** to learn about our professional pricing (In Canada please call 888-226-3042).



### Waterpik® Lunch & Learn Program

Water Pik, Inc. offers Lunch & Learn Programs for your office through our network of Professional Educators. To request a lunch and learn in your office please visit: [waterpik.com/lunch](http://waterpik.com/lunch)



Waterpik®



## FLUORIDE PRODUCTS

visit [professional.waterpik.com](http://professional.waterpik.com)  
or call **1.800.525.2020**

# waterpik® UltraThin & Clear™ Fluoride Varnish

### Glides on Smooth

- Contains 5% Sodium Fluoride
- Smooth, non-gritty texture
- Gluten Free
- Available in 4 flavors with Xylitol

#### Box of 30 single dose units 0.40mL and applicator brushes

20030674	Melon
20030675	Bubblegum
20030676	Strawberry
20030677	Mint

#### Box of 100 single dose units 0.40mL and applicator brushes

20030669	Melon
20030670	Bubblegum
20030671	Strawberry
20030672	Mint
20030673	Variety Pack (25 of each flavor)



## PROPHY PRODUCTS



### Waterpik® Soft Shine® Propy Paste Cleans and polishes in one step

No need to change pastes between restorations and natural teeth. Soft Shine® micro-fine white sapphire particles are proven to remove stains and fine scratches without damaging surfaces.

Safe and effective on restorations and natural teeth.

#### Jar of 100 cups

20008940 Ultra-Fine Grit Mint

## PROPHY PRODUCTS

### Waterpik® Prophy Pastes

Waterpik® prophylaxis pastes offer minimal splatter and great taste. Tough on stains and plaque, with minimal abrasion, this fluoride prophy paste makes cleaning and polishing effective for both patient and hygienist.

**Box of 200 single dose units (+1 prophy ring)**

- 20030445      Bubblegum Fine
- 20030446      Bubblegum Medium
- 20030455      Melon Fine
- 20030448      Melon Medium
- 20030449      Melon Coarse
- 20030450      Strawberry Medium
- 20030451      Strawberry Coarse
- 20030452      Cherry Coarse
- 20030453      Mint Extra Coarse
- 20030456      Mint Medium
- 20030454      Mint Coarse
- 20030447      Variety Pack Medium  
Includes 40 cups each:  
Bubblegum, Melon, Strawberry, Cherry, Mint



### Denco® Proflex Flexible Scaler

Designed for efficient scaling and stain removal, the Denco® proflex flexible scaler is made of stainless steel with microscopic scaling edges built into each strip. An indispensable tool for interproximal cleaning.

- 085000-000      3 foot roll



## WATERPIK® PROFESSIONAL DENTAL PRODUCTS SAMPLE PACK



To order, call  
**1-800-525-2020** or go to  
[waterpik.com/samplepack](http://waterpik.com/samplepack)

### Included in the Professional Dental Products Sample Pack:

**Waterpik® UltraThin & Clear Varnish .40ml single dose with applicator**  
**4 Flavors:**  
Melon  
Bubblegum  
Strawberry  
Mint

**Waterpik® Gel Grips Disposable Prophy Angles**  
1 - Soft Latex Free Dynamo Green Prophy cup (reverse helix design)  
1 - Natural Rubber Soft Blue Prophy cup (combined webbed and ribbed design)  
1 - Tapered Purple brush

**Waterpik® Prophy Paste single dose 1.23% Fluoride with Xylitol 7 Types:**  
Bubblegum Fine      Mint Medium  
Melon Medium      Mint Extra Coarse  
Strawberry Coarse      Cherry Coarse  
Mint Coarse

**Waterpik® Soft Shine 2gm. single dose**  
1 Ultra Fine Grit Mint

\*Not available in Canada

## DISPOSABLE PROPHY ANGLES

### waterpik Gel-Grip® Disposable Prophy Angle

#### Waterpik Gel-Grip® Prophy Angle with Latex-free Dynamo Green Cup

The cup features a unique reverse helix rib design to reduce splatter and efficiently remove stains. Soft cup makes contoured cleaning easy and maximizes patient comfort.

20018274      Soft Latex-Free Green prophy cup      200 angles

#### Waterpik Gel-Grip® Prophy Angles with Natural Rubber Cups

A combined webbed and ribbed cup design for stain removal and reduced splatter. Available in Soft Blue, for easy contour cleaning and patient comfort, and Firm Yellow for tougher stains.

20018272      Soft Blue prophy cup      200 angles

20018273      Firm Yellow prophy cup      200 angles

#### Waterpik Gel-Grip® Prophy Angle with Tapered Brush

Your best choice for removing heavy stains from occlusal surfaces and those hard-to-reach interproximal spaces.

20018276      Tapered Brush      100 angles



### waterpik® Densco® Disposable Prophy Angle

018125-100	Soft Latex-Free Green prophy cup	125 angles
018125-01M	Soft Latex-Free Green prophy cup	1,000 angles
018100-100	Soft Natural Rubber Blue prophy cup	125 angles
018100-01M	Soft Natural Rubber Blue prophy cup	1,000 angles
018150-100	Firm Natural Rubber Yellow prophy cup	125 angles
018150-01M	Firm Natural Rubber Yellow prophy cup	1,000 angles
20005474	Latex-Free cup/brush combination	125 angles
20005475	Latex-Free cup/brush combination	1,000 angles
018175-100	Tapered Brush	100 angles





## PROPHY CUPS




### Densco® Propphy Cups

Our complete line of prophy cups and brushes minimize splatter and fit most prophy angles. Densco® prophy cups are available in Soft Blue, Firm Yellow, Medium White, and Latex-free Dynamo Green, which features a reverse helix rib design for reduced splatter.



### Screw Shank Propphy Cups

- 085112-144 Soft Blue rubber with skirt, ribs only, 144 cups
- 085132-144 Soft Blue rubber with skirt, ribs and recessed webs, 144 cups
- 085132-01M Soft Blue rubber with skirt, ribs and recessed webs, 1000 cups
- 085133-144 Firm Yellow rubber with skirt, ribs and recessed webs, 144 cups
- 085133-01M Firm Yellow rubber with skirt, ribs and recessed webs, 1000 cups
- 085321-144 Medium White rubber with webs, 144 cups
- 085321-01M Medium White rubber with webs, 1000 cups



### Latex-free Dynamo Green Screw Shank Propphy Cups

- 085148-144 Soft Light Green rubber with reverse helix rib design, 144 cups
- 085148-01M Soft Light Green rubber with reverse helix rib design, 1000 cups



### Snap-On Propphy Cups

- 085032-144 Soft Blue rubber with ribs and recessed webs, 144 cups
- 085032-01M Soft Blue rubber with ribs and recessed webs, 1000 cups
- 085033-144 Firm Yellow rubber with ribs and recessed webs, 144 cups
- 085021-144 Medium White rubber with webs, 144 cups



### Latex-free Dynamo Green Snap-On Propphy Cups

- 085048-144 Soft Light Green rubber with reverse helix rib design, 144 cups



### Right Angle (Latch) Propphy Cups

- 085212-144 Soft Blue rubber with skirt, ribs only, 144 cups
- 085232-144 Soft Blue rubber with skirt, ribs and recessed webs, 144 cups
- 085232-01M Soft Blue rubber with skirt, ribs and recessed webs, 1000 cups
- 085233-144 Firm Yellow rubber with skirt, ribs and recessed webs, 144 cups
- 085233-01M Firm Yellow rubber with skirt, ribs and recessed webs, 1000 cups




### Densco® Propphy Brush

Densco® prophy brushes effectively remove stains and clean occlusal pits and fissures, around orthodontic appliances, and hard-to-reach interproximal spaces.



### Screw Shank Propphy Brush

- 085177-036 Tapered brush, 36 brushes



### Nylon Propphy Brush, Snap-on


- 085077-036 Tapered brush, 36 brushes

### Original Tofflemire® "Dead-Soft" Matrix Bands

The micro-thin "dead-soft" band is .001 inch gauge, making it more flexible than thicker bands. For use with posterior composites, it helps create tight interproximal contact. Allows creation of ideal restoration contour without breakage or band deformation.

**Packages of 36 bands**

- 34603-6 No.1 Adult Universal .001 in (.025mm)
- 34604-4 No.2 Adult MOD Wide .001 in (.025mm)
- 34605-2 No.13 Pedo Universal .001 in (.025mm)



<b>Thinnest</b>		
Micro-Thin	.001 in	.025 mm
<b>Thin</b>		
Ultra-Thin	.0015 in	.0375 mm
<b>Thickest</b>		
Regular	.002 in	.05 mm

### Original Tofflemire® Matrix Bands

This surgical grade stainless steel band is carefully stamped for consistently smooth (burr free) edges that keep your patients comfortable and reduce risk of gingival trauma. Available in a variety of sizes and contours, in both .002 inch and .0015 inch thicknesses.

**No.1 Adult Universal**

- 34577-0 .002 in (.05mm), 36 bands
- 062001-000 .002 in (.05mm), 144 bands
- 34578-9 .0015 in (.0375 mm), 36 bands
- 062002-000 .0015 in (.0375 mm), 144 bands
- 061032-000 .0015 in (.0375 mm), 1000 bands

**No.2 Adult MOD Wide**


- 34579-7 .002 in (.05mm), 36 bands
- 34580-0 .0015 in (.0375 mm), 36 bands

**No.3 Adult MOD Narrow**

- 34582-6 .0015 in (.0375 mm), 36 bands

**No.13 Pedo Universal**

- 34584-2 .002 in (.05mm), 36 bands
- 34585-0 .0015 in (.0375 mm), 36 bands
- 062010-000 .0015 in (.0375 mm), 144 bands



### Getz® Contour Matrix Bands "Dixieland Band"


The Getz® contoured matrix bands create anatomically contoured contacts with minimal burnishing for your best and tightest interproximal contacts. Made of high tensile strength stainless steel that will not tear under stress. Available in .002 inch regular and .0015 inch ultra-thin thickness.

**Packages of 50 bands**

- 061081-000 .002 in – Regular (.05mm)
- 061080-000 .0015 in – Ultra-thin (.0375mm)



## WATERPIK® WIZARD WEDGES®



Shaped for enhanced tooth separation

Exclusive multi-plane contour

### Wizard Wedges® Anatomical Wood Matrix Wedges

These anatomically shaped birch wood wedges employ an exclusive multi-plane contour that firmly engages adjoining teeth to hold the matrix in position. Available in four sizes.

**Bag of 100**

061315-000 Assortment of sizes

**Box of 400**

- 061306-000 Small (natural)
- 061307-000 Medium (yellow)
- 061308-000 Large (purple)
- 061309-000 Extra large (red)
- 061310-000 Assortment of sizes



### Wizard Wedges® Wood Matrix Wedges

These proven wedges ensure better access and stable holding of matrix systems during your critical restorative procedures. Specially crafted from birch wood for greater tooth separation and faster, better tooth position recovery upon removal. Available in four sizes.

**Package of 100**

061205-000 Assortment of sizes

**Box of 500**

- 061211-000 Small
- 061212-000 Large
- 061213-000 Slim Jim
- 061214-000 Wide Base
- 061215-000 Assortment of sizes



### Original Tofflemire® II Retainers

Secure locking mechanism and instant release enable the Original Tofflemire® retainers to tighten and release with precise, effortless action. The clean, smooth surface is free of burrs and sharp edges that can cause gingival trauma. Easy insertion and adjustment; handles any width band. Permits accurate trimming of wider band for special anatomical needs. Contra-Angle model can be used buccally and lingually. Fully autoclavable.

**Package of 1**

- 20016453 Universal
- 20016455 Contra-Angle



## IMPRESSION TRAYS

### Tra-Tens® Impression Trays

Proven over time to be a versatile design which provides extremely accurate impressions, Tra-Tens® impression trays have a unique, fully perforated pattern that enables excellent mechanical retention. Available in nine sizes.



**Packages of 12**

- 13451-000 No.1 Small Upper
- 13452-000 No.2 Medium Upper
- 13453-000 No.3 Large Upper
- 13454-000 No.4 Small Lower
- 13455-000 No.5 Medium Lower
- 13456-000 No.6 Large Lower
- 13457-000 No.7 Anterior
- 13458-000 No.8 Bilateral – Upper Left, Lower Right
- 13459-000 No.9 Bilateral – Upper Right, Lower Left

### Master Tray® Impression Trays

With a reinforced cross-section, Master Tray® disposable trays offer rigid metal tray-like performance for your most critical crown & bridge, and implant cases. Peripheral locking rim aids in retention. The Master Tray® is easy to adjust and trim. Available in nine sizes.

**Packages of 12**

- 011211-000 No.1 Large Upper
- 011212-000 No.2 Large Lower
- 011213-000 No.3 Medium Upper
- 011214-000 No.4 Medium Lower
- 011215-000 No.5 Small Upper
- 011216-000 No.6 Small Lower
- 011217-000 No.7 Partial – Upper Left, Lower Right
- 011218-000 No.8 Partial – Upper Right, Lower Left
- 011219-000 No.9 Anterior



### Sani-Trays® Dual Arch Impression Trays

Made with thin, strong mesh that will not interfere with occlusal bite. Available in Anterior, Posterior, Sideless Posterior, and Quadrant styles.

- Ⓐ 20000537 Anterior, box of 35
- Ⓑ 20000538 Posterior, box of 50
- Ⓒ 20000540 Sideless Posterior, box of 50
- Ⓓ 20000539 Quadrant, box of 35

### Sani-Trays® Impression Trays

An industry leader, Sani-Trays® impression trays produce precise alginate, crown and bridge, and edentulous impressions. Sani-Trays® feature extra thick sidewalls to minimize bending, even during complex procedures. Easy to adjust and trim. Available in 15 sizes.

**Packages of 12**

- 011511-012 No.1 Large Upper, Perforated
- 011512-012 No.2 Large Lower, Perforated
- 011513-012 No.3 Medium Upper, Perforated
- 011514-012 No.4 Medium Lower, Perforated
- 011515-012 No.5 Small Upper, Perforated
- 011516-012 No.6 Small Lower, Perforated
- 011517-012 No.7 Partial – Upper Left, Lower Right, Non-perforated
- 011518-012 No.8 Partial – Upper Right, Lower Left, Non-perforated
- 011519-012 No.9 Anterior, Non-perforated
- 011521-012 No.11 Large Upper, Non-perforated
- 011522-012 No.12 Large Lower, Non-perforated
- 011523-012 No.13 Medium Upper, Non-perforated
- 011524-012 No.14 Medium Lower, Non-perforated
- 011525-012 No.15 Small Upper, Non-perforated
- 011526-012 No.16 Small Lower, Non-perforated



## IMPRESSION MATERIALS



### Key-To® Alginate Impression Material

A versatile choice, Key-To® alginate comes in a variety of body and set types designed to produce consistent impressions. Quick mixing for smooth handling and fast cleanup. Withdraws easily from undercuts and bands without tearing or distortion. Each can includes powder scoop and water vial. Available in Regular and Fast Set. Spearmint flavored.

**Case of 10 one lb (454g) cans**

- 013111-000 Regular Body-Regular Set
- 013112-000 Heavy Body-Fast Set
- 013113-000 Regular Body-Fast Set

**25 lb drum (11.34kg)**

- 013126-000 Heavy Body-Fast Set

### Kalginate™ Alginate

Kalginate™ alginate is a heavier-bodied impression material specifically formulated to produce superior void-free impressions with minimal bubbles. Available in Normal and Fast-Setting formulas, each can includes powder scoop and water vial. Cinnamon flavored.

**Case of 10 one lb (454g) cans**

- 015021-000 Regular Set
- 015023-000 Fast Set



### Opotow® Standard Zinc Oxide Eugenol Impression Paste

Designed for use with full denture impressions in a custom tray or as a wash in a denture base. Blue color; Type I, Heavy-Bodied.

012200-000 Complete package contains "A" Hardener 2.8oz(80g), "B" Base 5.1oz(145g), mixing pad, and tube keys



### Krex® Zinc Oxide Eugenol Impression Corrective Paste

A zinc oxide eugenol impression paste that can be used for full denture impressions in a custom tray or as a wash. Type II, Medium-Bodied.

015100-000 Complete package contains "A" Hardener 2.8oz (80g), "B" Base 2.3oz (65g), mixing pad, and tube keys

015101-000 Economy package contains "A" Hardener 4.9oz (139g), "B" Base 4.1oz (116g), mixing pad, and tube keys

### Hold™ Tray Adhesive

Available in liquid or aerosol spray. Not recommended for use with elastomeric impression materials.

011460-000 .5oz (14ml) bottles, package of two  
 011462-000 4oz (118 ml) bottles, package of two  
 011461-000 2.1oz (80ml) aerosol can



## CEMENTS

### Neo-Temp™ Temporary Resin Cement

Neo-Temp™ resin cement is a non-eugenol temporary luting material based on a chemically cured urethane methacrylate polymer.

- 052200-000 Complete kit contains 10g cement ("A" syringe), 10g ("B" syringe), 15ml of dentin sealant, 1.5g releasing agent, brushes, spatulas, and mixing pad
- 052205-000 10g "A" syringe refill
- 052210-000 10g "B" syringe refill  
Not available in Canada



### Opotow® Temporary Cement

Designed for luting of cast crowns and bridges on a semi-permanent basis as well as for temporary cementation of provisional restorations where prolonged evaluation is required.

- 052117-000 Complete kit contains accelerator, base, mixing pad, and two tube keys  
Not available in Canada

### Opotow® Trial Cement

Opotow® trial cement is a non-hardening temporary cement designed for trial cementation of bridgework, splints, and porcelain-fused-to-metal crowns and bridges.

- 052120-000 Complete kit contains accelerator, base, mixing pad, and tube key  
Not available in Canada



## CROWN & BRIDGE

### Getz® Blue Core Crown Build-up Material

Getz® blue core crown build-up material is a small particle, radiopaque, light-cured composite core material that cures to a depth of 6mm in 30 seconds.

041418-000 Blue core crown build-up material, single 5g syringe



## DENTIN SEALANT



### Barrier Dentin Sealant

A universal protective cavity varnish that is also an ideal coronal sealant for endodontic procedures and an effective short-term protectant for exposed root surfaces.

042305-000 Complete kit contains 15ml of barrier dentin sealant and 50 disposable applicator pipettes  
 042315-000 Disposable applicator pipettes, 2<sup>7</sup>/<sub>8</sub> inch long, package of 50

## ACCESSORIES



### Getz® Composite Accessory

Replacement brushes for mixing and application of restorative materials.

Package of 100  
 041465-000 Disposable brushes



### Jiffy Tubes

Available in straight and curved models, Jiffy Tubes aid in the placement of cements, filling materials, and medicaments in difficult access areas.

Box of 30  
 051001-000 Straight tubes  
 051002-000 Curved tubes




## INSTRUMENTS

### Densco® Condensaire® Amalgam Condenser

The autoclavable Condensaire® amalgam condenser is a safe and effective way to condense amalgam automatically. Condensaire® is easy and safe to install and operate, fitting all standard four hole handpiece connectors. Requires no routine maintenance or special cleaning.

- 083117-000 Complete kit: Instrument with four hole connector (includes: Medium Round tip, Large Round tip, Round Rubber tip, and Tip holder)
- 083111-000 Instrument only with four hole connector



Holder sold separately

### Densco® Condensaire® Condenser Accessories

**Tips**

- 000809-000 Small Round (1.14mm)
- 000810-000 Medium Round (1.52mm)\*
- 000811-000 Large Round (2.16mm)\*
- 000812-000 Round Rubber (4.75mm)\*



- 000833-000 Rubber replacement tip, round
- 000845-000 Rubber replacement tip, pointed
- 000829-000 Tip holder\*

\*Included in Complete Kit (083117-000)

Water Pik, Inc. a subsidiary  
of Church & Dwight Co., Inc.  
1730 E. Prospect Road  
Fort Collins, CO 80553-0001

**waterpik.com**

Water Pik Professional Customer Service:

**1-800-525-2020**

In Canada, please call 1-888-226-3042



Sonicare®, Oral-B®, and Wizard Wedges®  
are not trademarks of Water Pik, Inc.

PN 20010232-CAT  
FN 20010232-CAT AU

© 2022 Water Pik, Inc.



**waterpik®**