

PROFESSIONAL ORAL HEALTH PRODUCT GUIDE

waterpik[®]

FLUORIDE PRODUCTS

visit professional.waterpik.com
or call 1.800.525.2020

waterpik® UltraThin & Clear™ Fluoride Varnish

Glides on Smooth

- Contains 5% Sodium Fluoride
- Smooth, non-gritty texture
- Gluten Free
- Available in 4 flavors with Xylitol

Box of 30 single dose units 0.40mL and applicator brushes

20030674	Melon
20030675	Bubblegum
20030676	Strawberry
20030677	Mint

Box of 100 single dose units 0.40mL and applicator brushes

20030669	Melon
20030670	Bubblegum
20030671	Strawberry
20030672	Mint
20030673	Variety Pack (25 of each flavor)



PROPHY PRODUCTS



Waterpik® Soft Shine® Propy Paste
Cleans and polishes in one step

No need to change pastes between restorations and natural teeth. Soft Shine® micro-fine white sapphire particles are proven to remove stains and fine scratches without damaging surfaces.
Safe and effective on restorations and natural teeth.

Jar of 100 cups
20008940

Ultra-Fine Grit Mint

PROPHY PRODUCTS

Waterpik® Prophecy Pastes

Waterpik® prophylaxis pastes offer minimal splatter and great taste. Tough on stains and plaque, with minimal abrasion, this fluoride prophylaxis paste makes cleaning and polishing effective for both patient and hygienist.

Box of 200 single dose units (+1 prophylaxis ring)

- 20030445 Bubblegum Fine
- 20030446 Bubblegum Medium
- 20030455 Melon Fine
- 20030448 Melon Medium
- 20030449 Melon Coarse
- 20030450 Strawberry Medium
- 20030451 Strawberry Coarse
- 20030452 Cherry Coarse
- 20030453 Mint Extra Coarse
- 20030456 Mint Medium
- 20030454 Mint Coarse
- 20030447 Variety Pack Medium
Includes 40 cups each:
Bubblegum, Melon, Strawberry, Cherry, Mint



Denco® Proflex Flexible Scaler

Designed for efficient scaling and stain removal, the Denco® proflex flexible scaler is made of stainless steel with microscopic scaling edges built into each strip. An indispensable tool for interproximal cleaning.

085000-000 3 foot roll



WATERPIK® PROFESSIONAL DENTAL PRODUCTS SAMPLE PACK



To order, call
1-800-525-2020 or go to
waterpik.com/samplepack

Included in the Professional Dental Products Sample Pack:

Waterpik® UltraThin & Clear Varnish .40ml single dose with applicator
4 Flavors:
Melon
Bubblegum
Strawberry
Mint

Waterpik® Gel Grips Disposable Propyl Angles
1 - Soft Latex Free Dynamo Green Propyl cup (reverse helix design)
1 - Natural Rubber Soft Blue Propyl cup (combined webbed and ribbed design)
1 - Tapered Purple brush

Waterpik® Propyl Paste single dose 1.23% Fluoride with Xylitol 7 Types:
Bubblegum Fine Mint Medium
Melon Medium Mint Extra Coarse
Strawberry Coarse Cherry Coarse
Mint Coarse

Waterpik® Soft Shine 2gm. single dose
1 Ultra Fine Grit Mint

*Not available in Canada

DISPOSABLE PROPHY ANGLES

waterpik® Gel-Grip® Disposable Prophy Angle

Waterpik Gel-Grip® Prophy Angle with Latex-free Dynamo Green Cup

The cup features a unique reverse helix rib design to reduce splatter and efficiently remove stains. Soft cup makes contoured cleaning easy and maximizes patient comfort.

20018274 Soft Latex-Free Green prophy cup 200 angles

Waterpik Gel-Grip® Prophy Angles with Natural Rubber Cups

A combined webbed and ribbed cup design for stain removal and reduced splatter. Available in Soft Blue, for easy contour cleaning and patient comfort, and Firm Yellow for tougher stains.

20018272 Soft Blue prophy cup 200 angles

20018273 Firm Yellow prophy cup 200 angles

Waterpik Gel-Grip® Prophy Angle with Tapered Brush

Your best choice for removing heavy stains from occlusal surfaces and those hard-to-reach interproximal spaces.

20018276 Tapered Brush 100 angles



waterpik® Densco® Disposable Prophy Angle

018125-100	Soft Latex-Free Green prophy cup	125 angles
018125-01M	Soft Latex-Free Green prophy cup	1,000 angles
018100-100	Soft Natural Rubber Blue prophy cup	125 angles
018100-01M	Soft Natural Rubber Blue prophy cup	1,000 angles
018150-100	Firm Natural Rubber Yellow prophy cup	125 angles
018150-01M	Firm Natural Rubber Yellow prophy cup	1,000 angles
20005474	Latex-Free cup/brush combination	125 angles
20005475	Latex-Free cup/brush combination	1,000 angles
018175-100	Tapered Brush	100 angles



PROPHY CUPS




Densco® Propphy Cups

Our complete line of prophy cups and brushes minimize splatter and fit most prophy angles. Densco® prophy cups are available in Soft Blue, Firm Yellow, Medium White, and Latex-free Dynamo Green, which features a reverse helix rib design for reduced splatter.



Screw Shank Propphy Cups

- 085112-144 Soft Blue rubber with skirt, ribs only, 144 cups
- 085132-144 Soft Blue rubber with skirt, ribs and recessed webs, 144 cups
- 085132-01M Soft Blue rubber with skirt, ribs and recessed webs, 1000 cups
- 085133-144 Firm Yellow rubber with skirt, ribs and recessed webs, 144 cups
- 085133-01M Firm Yellow rubber with skirt, ribs and recessed webs, 1000 cups
- 085321-144 Medium White rubber with webs, 144 cups
- 085321-01M Medium White rubber with webs, 1000 cups



Latex-free Dynamo Green Screw Shank Propphy Cups

- 085148-144 Soft Light Green rubber with reverse helix rib design, 144 cups
- 085148-01M Soft Light Green rubber with reverse helix rib design, 1000 cups




Snap-On Propphy Cups

- 085032-144 Soft Blue rubber with ribs and recessed webs, 144 cups
- 085032-01M Soft Blue rubber with ribs and recessed webs, 1000 cups
- 085033-144 Firm Yellow rubber with ribs and recessed webs, 144 cups
- 085021-144 Medium White rubber with webs, 144 cups



Latex-free Dynamo Green Snap-On Propphy Cups

- 085048-144 Soft Light Green rubber with reverse helix rib design, 144 cups



Right Angle (Latch) Propphy Cups

- 085212-144 Soft Blue rubber with skirt, ribs only, 144 cups
- 085232-144 Soft Blue rubber with skirt, ribs and recessed webs, 144 cups
- 085232-01M Soft Blue rubber with skirt, ribs and recessed webs, 1000 cups
- 085233-144 Firm Yellow rubber with skirt, ribs and recessed webs, 144 cups
- 085233-01M Firm Yellow rubber with skirt, ribs and recessed webs, 1000 cups




Densco® Propphy Brush

Densco® prophy brushes effectively remove stains and clean occlusal pits and fissures, around orthodontic appliances, and hard-to-reach interproximal spaces.



Screw Shank Propphy Brush

- 085177-036 Tapered brush, 36 brushes



Nylon Propphy Brush, Snap-on


- 085077-036 Tapered brush, 36 brushes

Original Tofflemire® "Dead-Soft" Matrix Bands

The micro-thin "dead-soft" band is .001 inch gauge, making it more flexible than thicker bands. For use with posterior composites, it helps create tight interproximal contact. Allows creation of ideal restoration contour without breakage or band deformation.

Packages of 36 bands


- 34603-6 No.1 Adult Universal .001 in (.025mm)
- 34604-4 No.2 Adult MOD Wide .001 in (.025mm)
- 34605-2 No.13 Pedo Universal .001 in (.025mm)



Thinnest Micro-Thin	.001 in	.025 mm
Thin Ultra-Thin	.0015 in	.0375 mm
Thickest Regular	.002 in	.05 mm

Original Tofflemire® Matrix Bands

This surgical grade stainless steel band is carefully stamped for consistently smooth (burr free) edges that keep your patients comfortable and reduce risk of gingival trauma. Available in a variety of sizes and contours, in both .002 inch and .0015 inch thicknesses.



No.1 Adult Universal

- 34577-0 .002 in (.05 mm), 36 bands
- 062001-000 .002 in (.05 mm), 144 bands
- 34578-9 .0015 in (.0375 mm), 36 bands
- 062002-000 .0015 in (.0375 mm), 144 bands
- 061032-000 .0015 in (.0375 mm), 1000 bands

No.2 Adult MOD Wide

- 34579-7 .002 in (.05 mm), 36 bands
- 34580-0 .0015 in (.0375 mm), 36 bands

No.3 Adult MOD Narrow

- 34582-6 .0015 in (.0375 mm), 36 bands

No.13 Pedo Universal

- 34584-2 .002 in (.05 mm), 36 bands
- 34585-0 .0015 in (.0375 mm), 36 bands
- 062010-000 .0015 in (.0375 mm), 144 bands

Getz® Contour Matrix Bands "Dixieland Band"


The Getz® contoured matrix bands create anatomically contoured contacts with minimal burnishing for your best and tightest interproximal contacts. Made of high tensile strength stainless steel that will not tear under stress. Available in .002 inch regular and .0015 inch ultra-thin thickness.

Packages of 50 bands

- 061081-000 .002 in – Regular (.05mm)
- 061080-000 .0015 in – Ultra-thin (.0375 mm)



WATERPIK® WIZARD WEDGE



Wizard Wedges® Anatomical Wood Matrix Wedges

These anatomically shaped birch wood wedges employ an exclusive multi-plane contour that firmly engages adjoining teeth to hold the matrix in position. Available in four sizes.

Bag of 100

061315-000 Assortment of sizes

Box of 400

061306-000 Small (natural)
 061307-000 Medium (yellow)
 061308-000 Large (purple)
 061309-000 Extra large (red)
 061310-000 Assortment of sizes



Wizard Wedges® Wood Matrix Wedges

These proven wedges ensure better access and stable holding of matrix systems during your critical restorative procedures. Specially crafted from birch wood for greater tooth separation and faster, better tooth position recovery upon removal. Available in four sizes.

Package of 100

061205-000 Assortment of sizes

Box of 500

061211-000 Small
 061212-000 Large
 061213-000 Slim Jim
 061214-000 Wide Base
 061215-000 Assortment of sizes



Original Tofflemire® II Retainers

Secure locking mechanism and instant release enable the Original Tofflemire® retainers to tighten and release with precise, effortless action. The clean, smooth surface is free of burrs and sharp edges that can cause gingival trauma. Easy insertion and adjustment; handles any width band. Permits accurate trimming of wider band for special anatomical needs. Contra-Angle model can be used buccally and lingually. Fully autoclavable.

Package of 1

20016453 Universal
 20016455 Contra-Angle



IMPRESSION TRAYS

Tra-Tens® Impression Trays

Proven over time to be a versatile design which provides extremely accurate impressions, Tra-Tens® impression trays have a unique, fully perforated pattern that enables excellent mechanical retention. Available in nine sizes.



Packages of 12

- 13451-000 No.1 Small Upper
- 13452-000 No.2 Medium Upper
- 13453-000 No.3 Large Upper
- 13454-000 No.4 Small Lower
- 13455-000 No.5 Medium Lower
- 13456-000 No.6 Large Lower
- 13457-000 No.7 Anterior
- 13458-000 No.8 Bilateral – Upper Left, Lower Right
- 13459-000 No.9 Bilateral – Upper Right, Lower Left

Master Tray® Impression Trays

With a reinforced cross-section, Master Tray® disposable trays offer rigid metal tray-like performance for your most critical crown & bridge, and implant cases. Peripheral locking rim aids in retention. The Master Tray® is easy to adjust and trim. Available in nine sizes.


Packages of 12

- 011211-000 No.1 Large Upper
- 011212-000 No.2 Large Lower
- 011213-000 No.3 Medium Upper
- 011214-000 No.4 Medium Lower
- 011215-000 No.5 Small Upper
- 011216-000 No.6 Small Lower
- 011217-000 No.7 Partial – Upper Left, Lower Right
- 011218-000 No.8 Partial – Upper Right, Lower Left
- 011219-000 No.9 Anterior



Sani-Trays® Dual Arch Impression Trays

Made with thin, strong mesh that will not interfere with occlusal bite. Available in Anterior, Posterior, Sideless Posterior, and Quadrant styles.



- (A) 20000537 Anterior, box of 35
- (B) 20000538 Posterior, box of 50
- (C) 20000540 Sideless Posterior, box of 50
- (D) 20000539 Quadrant, box of 35

Sani-Trays® Impression Trays

An industry leader, Sani-Trays® impression trays produce precise alginate, crown and bridge, and edentulous impressions. Sani-Trays® feature extra thick sidewalls to minimize bending, even during complex procedures. Easy to adjust and trim. Available in 15 sizes.

Packages of 12

011511-012	No.1 Large Upper, Perforated
011512-012	No.2 Large Lower, Perforated
011513-012	No.3 Medium Upper, Perforated
011514-012	No.4 Medium Lower, Perforated
011515-012	No.5 Small Upper, Perforated
011516-012	No.6 Small Lower, Perforated
011517-012	No.7 Partial – Upper Left, Lower Right, Non-perforated
011518-012	No.8 Partial – Upper Right, Lower Left, Non-perforated
011519-012	No.9 Anterior, Non-perforated
011521-012	No.11 Large Upper, Non-perforated
011522-012	No.12 Large Lower, Non-perforated
011523-012	No.13 Medium Upper, Non-perforated
011524-012	No.14 Medium Lower, Non-perforated
011525-012	No.15 Small Upper, Non-perforated
011526-012	No.16 Small Lower, Non-perforated



IMPRESSION MATERIALS



Key-To® Alginate Impression Material

A versatile choice, Key-To® alginate comes in a variety of body and set types designed to produce consistent impressions. Quick mixing for smooth handling and fast cleanup. Withdraws easily from undercuts and bands without tearing or distortion. Each can includes powder scoop and water vial. Available in Regular and Fast Set. Spearmint flavored.

Case of 10 one lb (454g) cans

013111-000	Regular Body-Regular Set
013112-000	Heavy Body-Fast Set
013113-000	Regular Body-Fast Set

25 lb drum (11.34kg)

013126-000	Heavy Body-Fast Set
------------	---------------------

Kalginate™ Alginate

Kalginate™ alginate is a heavier-bodied impression material specifically formulated to produce superior void-free impressions with minimal bubbles. Available in Normal and Fast-Setting formulas, each can includes powder scoop and water vial. Cinnamon flavored.

Case of 10 one lb (454g) cans

015021-000	Regular Set
015023-000	Fast Set



Opotow® Standard Zinc Oxide Eugenol Impression Paste

Designed for use with full denture impressions in a custom tray or as a wash in a denture base. Blue color; Type I, Heavy-Bodied.

012200-000 Complete package contains "A" Hardener 2.8oz(80g), "B" Base 5.1oz(145g), mixing pad, and tube keys



Krex® Zinc Oxide Eugenol Impression Corrective Paste

A zinc oxide eugenol impression paste that can be used for full denture impressions in a custom tray or as a wash. Type II, Medium-Bodied.

015100-000 Complete package contains "A" Hardener 2.8oz (80g), "B" Base 2.3oz (65g), mixing pad, and tube keys

015101-000 Economy package contains "A" Hardener 4.9oz (139g), "B" Base 4.1oz (116g), mixing pad, and tube keys

Hold™ Tray Adhesive

Available in liquid or aerosol spray. Not recommended for use with elastomeric impression materials.

011460-000 .5oz (14 ml) bottles, package of two
 011462-000 4oz (118ml) bottles, package of two
 011461-000 2.1oz (80ml) aerosol can



CEMENTS

Neo-Temp™ Temporary Resin Cement

Neo-Temp™ resin cement is a non-eugenol temporary luting material based on a chemically cured urethane methacrylate polymer.

- 052200-000 Complete kit contains 10g cement ("A" syringe), 10g ("B" syringe), 15ml of dentin sealant, 1.5g releasing agent, brushes, spatulas, and mixing pad
- 052205-000 10g "A" syringe refill
- 052210-000 10g "B" syringe refill
- Not available in Canada



Opatow® Temporary Cement

Designed for luting of cast crowns and bridges on a semi-permanent basis as well as for temporary cementation of provisional restorations where prolonged evaluation is required.

- 052117-000 Complete kit contains accelerator, base, mixing pad, and two tube keys
- Not available in Canada

Opatow® Trial Cement

Opatow® trial cement is a non-hardening temporary cement designed for trial cementation of bridgework, splints, and porcelain-fused-to-metal crowns and bridges.

- 052120-000 Complete kit contains accelerator, base, mixing pad, and tube key
- Not available in Canada



CROWN & BRIDGE

Getz® Blue Core Crown Build-up Material

Getz® blue core crown build-up material is a small particle, radiopaque, light-cured composite core material that cures to a depth of 6mm in 30 seconds.

041418-000 Blue core crown build-up material, single 5g syringe



DENTIN SEALANT



Barrier Dentin Sealant

A universal protective cavity varnish that is also an ideal coronal sealant for endodontic procedures and an effective short-term protectant for exposed root surfaces.

042305-000 Complete kit contains 15ml of barrier dentin sealant and 50 disposable applicator pipettes
 042315-000 Disposable applicator pipettes, 2 7/8 inch long, package of 50

ACCESSORIES



Getz® Composite Accessory

Replacement brushes for mixing and application of restorative materials.

Package of 100
 041465-000 Disposable brushes



Jiffy Tubes

Available in straight and curved models, Jiffy Tubes aid in the placement of cements, filling materials, and medicaments in difficult access areas.


Box of 30
 051001-000 Straight tubes
 051002-000 Curved tubes

INSTRUMENTS

Densco® Condensaire® Amalgam Condenser

The autoclavable Condensaire® amalgam condenser is a safe and effective way to condense amalgam automatically. Condensaire® is easy and safe to install and operate, fitting all standard four hole handpiece connectors. Requires no routine maintenance or special cleaning.

- 083117-000 Complete kit: Instrument with four hole connector (includes: Medium Round tip, Large Round tip, Round Rubber tip, and Tip holder)
- 083111-000 Instrument only with four hole connector



Holder sold separately

Densco® Condensaire® Condenser Accessories

Tips

- 000809-000 Small Round (1.14mm)
- 000810-000 Medium Round (1.52mm)*
- 000811-000 Large Round (2.16mm)*
- 000812-000 Round Rubber (4.75mm)*



- 000833-000 Rubber replacement tip, round
- 000845-000 Rubber replacement tip, pointed
- 000829-000 Tip holder*

*Included in Complete Kit (083117-000)

Water Pik, Inc., a subsidiary
of Church & Dwight Co., Inc.
1730 E. Prospect Road
Fort Collins, CO 80553-0001

waterpik.com

Water Pik Professional Customer Service:

1-800-525-2020

In Canada, please call 1-888-226-3042



Wizard Wedges® is not a trademark of Water Pik, Inc.

PN 20026877-CAT
FN 20026877-CAT AG

© 2022 Water Pik, Inc.
Made in the USA.



waterpik®